

Positope™

Catalog no. R900-50

General Information

Introduction

Positope™ is a recombinant protein specifically engineered to contain seven different tags for detection with seven different antibodies. Positope™ protein is intended for use as a positive control for antibody function in western blot experiments.

Shipping/Storage

The Positope™ protein is shipped at room temperature. Upon receipt, please store at +4°C. The protein is stable for 6 months at +4°C.

Contents

Positope™ recombinant protein is supplied in a total volume of 200 µl at a concentration of 25 ng/µl in reducing SDS-PAGE sample buffer (63 mM Tris-HCl, 10% glycerol, 2% SDS, 0.0025% Bromophenol Blue, 50 mM β-mercaptoethanol).

Product Qualification

Each lot of Positope™ must be purified to > 95% purity from an *E. coli* lysate. In addition, each primary antibody listed below will detect 250 ng of Positope™ on a western blot using HRP-conjugated secondary antibody and chemiluminescence reagents. A strong signal is observed after a 1 minute exposure.

Description

Positope™ is a 53 kDa recombinant protein created using the gene for green fluorescent protein (GFP) and engineered to contain six additional epitopes as shown in the figure. It is expressed from pBAD/Thio-TOPO® in TOP10 cells using arabinose as an inducer. The table describes the actual epitope if known. Amino acid sequence is available by contacting Technical Service (page 3).

Note: HRP conjugates of Invitrogen antibodies, i.e. Anti-myc-HRP, Anti-V5-HRP, Anti-Xpress™-HRP, Anti-HisG-HRP, and Anti-His(C-term)-HRP, also detect Positope™.



Tag	Epitope	Reference
Thioredoxin	Unknown	(Dickason <i>et al.</i> , 1995)
HisG	-HHHHHHG-	(Robson <i>et al.</i> , 1995)
Xpress™	-DLYDDDDK-	(Kroll <i>et al.</i> , 1993)
GFP	Many	Invitrogen
c-myc	-EQKLISEEDL-	(Evan <i>et al.</i> , 1985)
V5	-GKPIPPLLGLDST-	(Southern <i>et al.</i> , 1991)
His(C-term)	-HHHHHH-COOH	(Lindner <i>et al.</i> , 1997)

Using Positope™



Please note that the Positope™ control protein is provided in reducing SDS-PAGE sample buffer and is **not** suitable for use in immunoassays (i.e. ELISA).

Western Analysis

In experiments performed at Invitrogen, we tested various amounts of Positope™ to obtain a strong signal in a western blot. Results are shown in the table below. Conditions were:

Primary antibody: 2 µl primary antibody in 20 ml buffer (1:10,000)

HRP-conjugated secondary antibody (goat anti-mouse-HRP or goat anti-rabbit-HRP):
1:10,000 dilution

Detection: Chemiluminescence reagents (ECL; Amersham)

Exposure: 1 minute

Primary Antibody	Amount of Positope™
Anti-Thio	100 ng
Anti-HisG	100 ng
Anti-Xpress™	250 ng
Anti-GFP (polyclonal)	100 ng
Anti-Myc	250 ng
Anti-V5	100 ng
Anti-His(C-term)	250 ng

Recommended Use

For Western blot analysis, we recommend using 250 ng (10 µl) of Positope™ per lane. Lower amounts of Positope™ may be used, but the strength of the signal may be affected. Positope™ recombinant protein is supplied in reducing SDS-PAGE sample buffer. We recommend that you perform the following steps before use:

1. Transfer the appropriate amount of Positope™ that you will load onto the gel (e.g. 10 µl) to a separate microcentrifuge tube.
2. Treat the Positope™ in the same manner as your other samples prior to loading your gel. We generally boil the sample for five minutes.



Note

Increasing the amount of Positope™ loaded, using too much antibody, or increasing the detection/exposure time may lead to the detection of proteolytic breakdown products of Positope™, especially if you are using the Anti-V5 Antibody.

For More Information

For more information on your specific antibody or for procedures to perform western blot analysis, please refer to the particular manual for your antibody.

Technical Service

World Wide Web



Visit the [Invitrogen Web Resource](#) using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your web browser (Netscape 3.0 or newer), then enter the following location (or URL):

<http://www.invitrogen.com>

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

Contact us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our web page (www.invitrogen.com).

United States Headquarters:

Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail:
tech_service@invitrogen.com

Japanese Headquarters

Invitrogen Japan K.K.
Nihonbashi Hama-Cho Park Bldg. 4F
2-35-4, Hama-Cho, Nihonbashi
Tel: 81 3 3663 7972
Fax: 81 3 3663 8242
E-mail: jpinfo@invitrogen.com

European Headquarters:

Invitrogen Ltd
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Tel (Free Phone Orders): 0800 269 210
Tel (General Enquiries): 0800 5345 5345
Fax: +44 (0) 141 814 6287
E-mail: eurotech@invitrogen.com

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

References

- Dickason, R. R., Edward, R. A., Bryan, J., and Huston, D. P. (1995). Versatile *E. coli* Thioredoxin Specific Monoclonal Antibodies Afford Convenient Analysis and Purification of Prokaryotic Expressed Soluble Fusion Protein. *J. Immunol. Methods* 185, 237-244.
- Evan, G. I., Lewis, G. K., Ramsay, G., and Bishop, V. M. (1985). Isolation of Monoclonal Antibodies Specific for *c-myc* Proto-oncogene Product. *Mol. Cell. Biol.* 5, 3610-3616.
- Kroll, D. J., Abdel-Hafiz, H. A.-M., Marcell, T., Simpson, S., Chen, C.-Y., Gutierrez-Hartmann, A., Lustbader, J. W., and Hoeffler, J. (1993). A Multifunctional Prokaryotic Protein Expression System: Overproduction, Affinity Purification, and Selective Detection. *DNA and Cell Biology* 12, 441-443.
- Lindner, P., Bauer, K., Krebber, A., Nieba, L., Kremmer, E., Krebber, C., Honegger, A., Klinger, B., Mocikat, R., and Pluckthun, A. (1997). Specific Detection of His-tagged Proteins With Recombinant Anti-His Tag scFv-Phosphatase or scFv-Phage Fusions. *BioTechniques* 22, 140-149.
- Robson, K. J. H., Frevert, U., Reckmann, I., Cowan, G., Beier, J., Scragg, I. G., Takehara, K., Bishop, D. H. L., Pradel, G., Sinden, R., Saccheo, S., Muller, H.-M., and Crisanti, A. (1995). Thrombospondin-Related Adhesive Protein (TRAP) of *Plasmodium falciparum*: Expression During Sporozite Ontogeny and Binding to Human Hepatocytes. *EMBO J.* 14, 3883-3894.
- Southern, J. A., Young, D. F., Heaney, F., Baumgartner, W., and Randall, R. E. (1991). Identification of an Epitope on the P and V Proteins of Simian Virus 5 That Distinguishes Between Two Isolates with Different Biological Characteristics. *J. Gen. Virol.* 72, 1551-1557.

©1999-2001 Invitrogen Corporation. All rights reserved.