


Method

Method ID	pH_calibration_7_4_10	Name	pH_calibration_7_4_10
Internal ID	M11	Version	14
Valid	Yes	Release state	Released
User defined label 1			
Can run on blocked instr.	No		

Comment



Export templates used

pH_cal_raw_resource_tabl CSV Simple
e_data

Sample (Calib) (C1)

Sensor

Sensor type	pH
Sensor	DG115-SC
Unit	pH

Calibration

Action	Calibration
Automatic buffer recognition	No
pH buffer list	FISHER (Ref. 25°C)
Calibration	Segmented
Number of buffers	3
Buffer 1	7.00
Buffer 2	4.00
Buffer 3	10.00

Temperature

Temperature [°C]	25.0
------------------	------

Titration stand (Titrationstand1)

Titration stand

Type	Rondolino TTL
Titration stand	Rondolino TTL 1

Stir (Stir1)

Stir

Speed [%]	30
Duration [s]	60
Condition	No

Measure (normal) [1] (Measurenormal1)

Sensor

Type	pH
Sensor	DG115-SC
Unit	mV

Temperature acquisition

Temperature acquisition	Yes
Temperature sensor	DT1000
Temperature unit	°C

Stir

Speed [%]	30
-----------	----

Acquisition of measured values

Acquisition	Equilibrium controlled
dE [mV]	0.5
dt [s]	1
t(min) [s]	30
t(max) [s]	60
Mean value	No

Condition

Condition	No
-----------	----

Calibration (Calibration1)

Calibration

Sensor	DG115-SC
Limits	No
Condition	No

Record (Record2)

Report

Report template	Calibration & titration report
Print	No
Condition	No

Signing History

Time	Full name	User login name	Role name	Operation	Comment